

产品名称: CBX4 Rabbit Monoclonal Antibody
产品货号: AMRe02830



产品概述 (Summary)

产品名称 (Production Name)	CBX4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CBX4
别名 (Alternative Names)	CBX4; hPc2; NBP16; PC2
基因 ID (Gene ID)	8535
蛋白 ID (SwissProt ID)	O00257.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 61 kDa; Observed MW: 61 kDa

研究背景 (Background)

E3 SUMO-protein ligase which facilitates SUMO1 conjugation by UBE2I. Component of the Polycomb group (PcG)

产品名称: CBX4 Rabbit Monoclonal Antibody
产品货号: AMRe02830




multiprotein PRC1 complex, a complex required to maintain the transcriptionally repressive state of many genes, including Hox genes, throughout development.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of CBX4 in C6, Raw264.7 lysates using CBX4 antibody.



Western blot analysis of CBX4 in HeLa, A549, HL-60, U251 lysates using CBX4 antibody.

注意事项 (Note)

For research use only .