

产品名称: Glucocorticoid Receptor Rabbit Monoclonal Antibody
产品货号: AMRe02833



产品概述 (Summary)

产品名称 (Production Name)	Glucocorticoid Receptor Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat, Hamster

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NR3C1
别名 (Alternative Names)	NR3C1; GRL; Glucocorticoid receptor; GR; Nuclear receptor subfamily 3 group C member 1
基因 ID (Gene ID)	2908
蛋白 ID (SwissProt ID)	P04150.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 86 kDa; Observed MW: 94 kDa

产品名称: Glucocorticoid Receptor Rabbit Monoclonal Antibody
产品货号: AMRe02833




研究背景 (Background)

The protein encoded by this gene is a receptor for glucocorticoids and can act as both a transcription factor and a regulator of other transcription factors. The encoded protein can bind DNA as a homodimer or as a heterodimer with another protein such as the retinoid X receptor. This protein can also be found in heteromeric cytoplasmic complexes along with heat shock factors and immunophilins.


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of Glucocorticoid Receptor in 3T3, Raw264.7 lysates using Glucocorticoid Receptor antibody.



Western blot analysis of Glucocorticoid Receptor in CHO-K1, C6, rat Brain, Ramos lysates using Glucocorticoid Receptor antibody.

注意事项 (Note)

For research use only .