

**产品名称: Phospho-Ezrin/Radixin/Moesin
(Thr567/Thr564/Thr558) Rabbit Monoclonal Antibody**
产品货号: AMRe02850



产品概述 (Summary)

产品名称 (Production Name)	Phospho-Ezrin/Radixin/Moesin (Thr567/Thr564/Thr558) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat,Hamster

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EZR
别名 (Alternative Names)	EZR; VIL2; Ezrin; Cytovillin; Villin-2; p81
基因 ID (Gene ID)	7430
蛋白 ID (SwissProt ID)	P15311.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 69 kDa; Observed MW: 75,80 kDa

**产品名称: Phospho-Ezrin/Radixin/Moesin
(Thr567/Thr564/Thr558) Rabbit Monoclonal Antibody**
产品货号: AMRe02850




研究背景 (Background)

Probably involved in connections of major cytoskeletal structures to the plasma membrane. In epithelial cells, required for the formation of microvilli and membrane ruffles on the apical pole. Along with PLEKHG6, required for normal macropinocytosis.


研究领域 (Research Area)

Cell Biology

图片 (Image Data)




Western blot analysis of Phospho-Ezrin (Thr567)/ Radixin (Thr564)/ Moesin (Thr558) in 3T3, 3T3+PDGF lysates using Phospho-Ezrin/Radixin/Moesin (Thr567/Thr564/Thr558) antibody.



Western blot analysis of Ezrin (Phospho-Thr567)/ Radixin (Phospho-Thr564)/ Moesin (Phospho-Thr558) in HeLa, CHO-K1, C6 lysates using Ezrin (Phospho-Thr567)/ Radixin (Phospho-Thr564)/ Moesin (Phospho-Thr558) antibody.

**产品名称: Phospho-Ezrin/Radixin/Moesin
(Thr567/Thr564/Thr558) Rabbit Monoclonal Antibody**
产品货号: AMRe02850



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Ezrin (Phospho-Thr567)/ Radixin (Phospho-Thr564)/ Moesin (Phospho-Thr558) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .