

产品名称: 14-3-3 alpha/beta Rabbit Monoclonal Antibody
产品货号: AMRe02880



产品概述 (Summary)

产品名称 (Production Name)	14-3-3 alpha/beta Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	YWHAB
别名 (Alternative Names)	YWHAB; 14-3-3 protein beta/alpha; Protein 1054; Protein kinase C inhibitor protein 1; KCIP-1
基因 ID (Gene ID)	7529
蛋白 ID (SwissProt ID)	P31946.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 28 kDa; Observed MW: 28 kDa

研究背景 (Background)

产品名称: 14-3-3 alpha/beta Rabbit Monoclonal Antibody
产品货号: AMRe02880




Adapter protein implicated in the regulation of a large spectrum of both general and specialized signaling pathways. Binds to a large number of partners, usually by recognition of a phosphoserine or phosphothreonine motif. Binding generally results in the modulation of the activity of the binding partner. Negative regulator of osteogenesis.


研究领域 (Research Area)

Neuroscience


图片 (Image Data)



Western blot analysis of 14-3-3 alpha/beta in rat brain, rat lung lysates using 14-3-3 alpha/beta antibody.



Western blot analysis of 14-3-3 alpha/beta in Hela, A549, MCF-7, U251, C2C12 lysates using 14-3-3 alpha/beta antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using 14-3-3 alpha/beta antibody. High-pressure and

产品名称: 14-3-3 alpha/beta Rabbit Monoclonal Antibody
产品货号: AMRe02880



temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .