

产品名称: ADK Rabbit Monoclonal Antibody
产品货号: AMRe02890



产品概述 (Summary)

产品名称 (Production Name)	ADK Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ADK
别名 (Alternative Names)	ADK; Adenosine kinase; AK; Adenosine 5'-phosphotransferase
基因 ID (Gene ID)	132
蛋白 ID (SwissProt ID)	P55263.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 41 kDa; Observed MW: 41 kDa

研究背景 (Background)

ATP dependent phosphorylation of adenosine and other related nucleoside analogs to monophosphate derivatives. Serves

产品名称: ADK Rabbit Monoclonal Antibody
产品货号: AMRe02890




as a potential regulator of concentrations of extracellular adenosine and intracellular adenine nucleotides.


研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Western blot analysis of ADK in PC-12, Raw264.7 lysates using ADK antibody.



Western blot analysis of ADK in Hela, A549, HL-60, U251, U87-MG lysates using ADK antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using ADK antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: ADK Rabbit Monoclonal Antibody
产品货号: AMRe02890



注意事项 (Note)

For research use only .