

产品名称: AMSH Rabbit Monoclonal Antibody
产品货号: AMRe02897



产品概述 (Summary)

产品名称 (Production Name)	AMSH Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	STAMBP
别名 (Alternative Names)	AMSH; MICCAP
基因 ID (Gene ID)	10617
蛋白 ID (SwissProt ID)	O95630.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 48 kDa; Observed MW: 48 kDa

研究背景 (Background)

Cytokine-mediated signal transduction in the JAK-STAT cascade requires the involvement of adaptor molecules. One such

产品名称: AMSH Rabbit Monoclonal Antibody
产品货号: AMRe02897




signal-transducing adaptor molecule contains an SH3 domain that is required for induction of MYC and cell growth. The protein encoded by this gene binds to the SH3 domain of the signal-transducing adaptor molecule, and plays a critical role in cytokine-mediated signaling for MYC induction and cell cycle progression. Multiple alternatively spliced transcript variants encoding the same protein isoform have been found for this gene.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of AMSH in A549, HL-60, U2OS lysates using AMSH antibody.

注意事项 (Note)

For research use only .