

产品名称: BCKDK Rabbit Monoclonal Antibody
产品货号: AMRe02910



产品概述 (Summary)

产品名称 (Production Name)	BCKDK Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BCKDK
别名 (Alternative Names)	BDK; BCKDKD
基因 ID (Gene ID)	10295
蛋白 ID (SwissProt ID)	O14874.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 46 kDa; Observed MW: 46 kDa

研究背景 (Background)

Catalyzes the phosphorylation and inactivation of the branched-chain alpha-ketoacid dehydrogenase complex, the key

产品名称: BCKDK Rabbit Monoclonal Antibody
产品货号: AMRe02910




regulatory enzyme of the valine, leucine and isoleucine catabolic pathways.


研究领域 (Research Area)

Tags & Cell Markers; Subcellular Markers; Organelles; Mitochondria; Signal Transduction; Metabolism; Amino Acids; Mitochondrial; Metabolism; Pathways and Processes; Mitochondrial Metabolism; Mitochondrial markers; Amino acid metabolism

图片 (Image Data)



Western blot analysis of BCKDK in PC-12, Raw264.7 lysates using BCKDK antibody.



Western blot analysis of BCKDK in HeLa, A549, HL-60 lysates using BCKDK antibody.

注意事项 (Note)

For research use only .