

产品名称: Integrin alpha L Rabbit Monoclonal Antibody
产品货号: AMRe02932



产品概述 (Summary)

产品名称 (Production Name)	Integrin alpha L Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Itgal
别名 (Alternative Names)	Cd11a; LFA-1; Ly-15; Ly-21; (p180); LFA-1A
基因 ID (Gene ID)	16408.0
蛋白 ID (SwissProt ID)	P24063.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 128 kDa; Observed MW: 160 kDa

研究背景 (Background)

Integrin ITGAL/ITGB2 is a receptor for ICAM1, ICAM2, ICAM3 and ICAM4.


产品名称: Integrin alpha L Rabbit Monoclonal Antibody
产品货号: AMRe02932




研究领域 (Research Area)

Tags & Cell Markers

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded mouse lung using Integrin alpha L antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of CD11a in A549, HL-60 lysates using CD11a antibody.

注意事项 (Note)

For research use only .