

产品名称: CrkL Rabbit Monoclonal Antibody
产品货号: AMRe02952



产品概述 (Summary)

产品名称 (Production Name)	CrkL Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CRKL
别名 (Alternative Names)	Crk L; Crk like protein; Crkl; Crkol
基因 ID (Gene ID)	1399
蛋白 ID (SwissProt ID)	P46109.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 34 kDa; Observed MW: 39 kDa

研究背景 (Background)

May mediate the transduction of intracellular signals.


产品名称: CrkL Rabbit Monoclonal Antibody
产品货号: AMRe02952




研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Immunocytochemistry analysis of CrkL (green) in K562 using CrkL antibody, and DAPI(blue).



Western blot analysis of CrkL in Jurkat lysates using CrkL antibody.

注意事项 (Note)

For research use only .