

产品名称: CTBP1 Rabbit Monoclonal Antibody
产品货号: AMRe02956



产品概述 (Summary)

产品名称 (Production Name)	CTBP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CTBP1
别名 (Alternative Names)	CTBP1; CTBP; C-terminal-binding protein 1; CtBP1
基因 ID (Gene ID)	1487
蛋白 ID (SwissProt ID)	Q13363.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 48 kDa; Observed MW: 48 kDa

研究背景 (Background)

Involved in controlling the equilibrium between tubular and stacked structures in the Golgi complex. Functions in brown

产品名称: CTBP1 Rabbit Monoclonal Antibody
产品货号: AMRe02956




adipose tissue (BAT) differentiation. Corepressor targeting diverse transcription regulators such as GLIS2. Has dehydrogenase activity.


研究领域 (Research Area)

Cell Biology


图片 (Image Data)



Immunocytochemistry analysis of CTBP1 (green) in HEPG2 using CTBP1 antibody, and DAPI (blue).



Western blot analysis of CtBP1 in A549, MCF-7, U251 lysates using CtBP1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using CtBP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: CTBP1 Rabbit Monoclonal Antibody
产品货号: AMRe02956



注意事项 (Note)

For research use only .