

产品名称: CUG BP1 Rabbit Monoclonal Antibody
产品货号: AMRe02958



产品概述 (Summary)

产品名称 (Production Name)	CUG BP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CELF1 CELF1; BRUNOL2; CUGBP; CUGBP1; NAB50; CUGBP Elav-like family member 1;
别名 (Alternative Names)	CELF-1; 50 kDa nuclear polyadenylated RNA-binding protein; Bruno-like protein 2; CUG triplet repeat RNA-binding protein 1; CUG-BP1; CUG-BP- and ETR-3-like factor 1; Dead
基因 ID (Gene ID)	10658
蛋白 ID (SwissProt ID)	Q92879.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 52 kDa; Observed MW: 52 kDa

产品名称: CUG BP1 Rabbit Monoclonal Antibody
产品货号: AMRe02958




研究背景 (Background)


RNA-binding protein implicated in the regulation of several post-transcriptional events. Involved in pre-mRNA alternative splicing, mRNA translation and stability. Mediates exon inclusion and/or exclusion in pre-mRNA that are subject to tissue-specific and developmentally regulated alternative splicing.

研究领域 (Research Area)

图片 (Image Data)




Immunocytochemistry analysis of CUG BP1 (green) in HT-1080 using CUG BP1 antibody, and DAPI (blue).




Western blot analysis of CUGBP1 in rat brain lysates using CUGBP1 antibody.

产品名称: CUG BP1 Rabbit Monoclonal Antibody
产品货号: AMRe02958



Western blot analysis of CUGBP1 in HeLa, Ramos lysates using CUGBP1 antibody



Immunohistochemistry analysis of paraffin-embedded Human tonsil using CUGBP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .