

产品名称: Cytokeratin 9 Rabbit Monoclonal Antibody
产品货号: AMRe02964



产品概述 (Summary)

产品名称 (Production Name)	Cytokeratin 9 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KRT9
别名 (Alternative Names)	K9; CK-9; EPPK
基因 ID (Gene ID)	3857
蛋白 ID (SwissProt ID)	P35527.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 62 kDa; Observed MW: 62 kDa

研究背景 (Background)

May serve an important special function either in the mature palmar and plantar skin tissue or in the morphogenetic

产品名称: Cytokeratin 9 Rabbit Monoclonal Antibody
产品货号: AMRe02964




program of the formation of these tissues. Plays a role in keratin filament assembly.


研究领域 (Research Area)

Cell Biology


图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded mouse skin using Cytokeratin 9 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Cytokeratin 9 in mouse brain lysates using Cytokeratin 9 antibody.



Western blot analysis of Cytokeratin 9 in C6, Ramos lysates using Cytokeratin 9 antibody

产品名称: Cytokeratin 9 Rabbit Monoclonal Antibody
产品货号: AMRe02964



注意事项 (Note)

For research use only .