

产品名称: Daxx Rabbit Monoclonal Antibody
产品货号: AMRe02966



产品概述 (Summary)

产品名称 (Production Name)	Daxx Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DAXX
别名 (Alternative Names)	DAXX; BING2; DAP6; Death domain-associated protein 6; Daxx; hDaxx; ETS1-associated protein 1; EAP1; Fas death domain-associated protein
基因 ID (Gene ID)	1616
蛋白 ID (SwissProt ID)	Q9UER7.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 81 kDa; Observed MW: 110 kDa

研究背景 (Background)

产品名称: Daxx Rabbit Monoclonal Antibody
产品货号: AMRe02966




Acts as an adapter protein in a MDM2-DAXX-USP7 complex by regulating the RING-finger E3 ligase MDM2 ubiquitination activity. Under non-stress condition, in association with the deubiquitinating USP7, prevents MDM2 self-ubiquitination and enhances the intrinsic E3 ligase activity of MDM2 towards TP53, thereby promoting TP53 ubiquitination and subsequent proteasomal degradation. Upon DNA damage, its association with MDM2 and USP7 is disrupted.


研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Immunocytochemistry analysis of Daxx (green) in K562 using Daxx antibody, and DAPI(blue).



Western blot analysis of Daxx in HepG2, U251, HL-60 lysates using Daxx antibody.

注意事项 (Note)

For research use only .