

产品名称: DNAJC12 Rabbit Monoclonal Antibody
产品货号: AMRe02973



产品概述 (Summary)

产品名称 (Production Name)	DNAJC12 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DNAJC12
别名 (Alternative Names)	JDP1; HPANBH4
基因 ID (Gene ID)	56521
蛋白 ID (SwissProt ID)	Q9UKB3.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 23 kDa; Observed MW: 23 kDa

研究背景 (Background)

This gene encodes a member of a subclass of the HSP40/DnaJ protein family. Members of this family of proteins are

产品名称: DNAJC12 Rabbit Monoclonal Antibody
产品货号: AMRe02973




associated with complex assembly, protein folding, and export. Two transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of DNAJC12 in HeLa, A549, HL-60 lysates using DNAJC12 antibody.

注意事项 (Note)

For research use only .