

产品名称: Elongation Factor 1A1 Rabbit Monoclonal Antibody
产品货号: AMRe02977



产品概述 (Summary)

产品名称 (Production Name)	Elongation Factor 1A1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EEF1A1
别名 (Alternative Names)	EEF1A1; CCS-3; CCS3; EEF-1; EEF1A; EF-Tu; EF1A; GRAF-1EF; LENG7; PTI1; eEF1A-1
基因 ID (Gene ID)	1915
蛋白 ID (SwissProt ID)	P68104.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 50 kDa; Observed MW: 50 kDa

产品名称: Elongation Factor 1A1 Rabbit Monoclonal Antibody
产品货号: AMRe02977




研究背景 (Background)

This protein promotes the GTP-dependent binding of aminoacyl-tRNA to the A-site of ribosomes during protein biosynthesis.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)




Immunocytochemistry analysis of Elongation Factor 1A1 (green) in HeLa using Elongation Factor 1A1 antibody, and DAPI(blue).



Western blot analysis of eEF1A1/EFTu in U2OS, C6 lysates using eEF1A1/EFTu antibody.

产品名称: Elongation Factor 1A1 Rabbit Monoclonal Antibody
产品货号: AMRe02977



Immunohistochemistry analysis of paraffin-embedded Human tonsil using eEF1A1/EFTu antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .