

产品名称: Entactin Rabbit Monoclonal Antibody
产品货号: AMRe02984



产品概述 (Summary)

产品名称 (Production Name)	Entactin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NID1
别名 (Alternative Names)	Enactin; Ent; Entactin 1; Entactin; NID 1; NID; Nid1; Nidogen 1; Nidogen
基因 ID (Gene ID)	4811
蛋白 ID (SwissProt ID)	P14543.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 136 kDa; Observed MW: 136 kDa

研究背景 (Background)

Sulfated glycoprotein widely distributed in basement membranes and tightly associated with laminin. Also binds to

产品名称: Entactin Rabbit Monoclonal Antibody
产品货号: AMRe02984




collagen IV and perlecan. It probably has a role in cell-extracellular matrix interactions.


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of Entactin/NID in Hela, A549, HL-60, U251, U87-MG lysates using Entactin antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Entactin/NID antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .