

产品名称: EPS15R Rabbit Monoclonal Antibody
产品货号: AMRe02987



产品概述 (Summary)

产品名称 (Production Name)	EPS15R Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EPS15L1
别名 (Alternative Names)	EPS15R
基因 ID (Gene ID)	58513
蛋白 ID (SwissProt ID)	Q9UBC2.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 94 kDa; Observed MW: 130 kDa

研究背景 (Background)

Seems to be a constitutive component of clathrin-coated pits that is required for receptor-mediated endocytosis. Involved

产品名称: EPS15R Rabbit Monoclonal Antibody
产品货号: AMRe02987




in endocytosis of integrin beta-1 (ITGB1) and transferrin receptor (TFR); internalization of ITGB1 as DAB2-dependent cargo but not TFR seems to require association with DAB2.


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Immunocytochemistry analysis of EPS15R (green) in HeLa using EPS15R antibody, and DAPI (blue).



Western blot analysis of EPS15R in A549, HL-60 lysates using EPS15R antibody.

注意事项 (Note)

For research use only .