

产品名称: ErbB 4 Rabbit Monoclonal Antibody
产品货号: AMRe02989



产品概述 (Summary)

产品名称 (Production Name)	ErbB 4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ERBB4
别名 (Alternative Names)	ERBB4; HER4; Receptor tyrosine-protein kinase erbB-4; Proto-oncogene-like protein c-ErbB-4; Tyrosine kinase-type cell surface receptor HER4; p180erbB4
基因 ID (Gene ID)	2066
蛋白 ID (SwissProt ID)	Q15303.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 147 kDa; Observed MW: 180 kDa

研究背景 (Background)

产品名称: ErbB 4 Rabbit Monoclonal Antibody
产品货号: AMRe02989




The HER4/ERBB4 gene is a member of the type I receptor tyrosine kinase subfamily that includes EGFR (MIM 131550), ERBB2 (MIM 164870), and ERBB3 (MIM 190151). It encodes a receptor for NDF/heregulin (MIM 142445).

研究领域 (Research Area)

Cancer

图片 (Image Data)



Western blot analysis of ErbB 4 in HeLa, A549, HL-60, U2OS, C6 lysates using ErbB 4 antibody.

注意事项 (Note)

For research use only .