

产品名称: Ezrin Rabbit Monoclonal Antibody
产品货号: AMRe02993



产品概述 (Summary)

产品名称 (Production Name)	Ezrin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EZR
别名 (Alternative Names)	EZR; VIL2; Ezrin; Cytovillin; Villin-2; p81
基因 ID (Gene ID)	7430
蛋白 ID (SwissProt ID)	P15311.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 69 kDa; Observed MW: 80 kDa

研究背景 (Background)

Probably involved in connections of major cytoskeletal structures to the plasma membrane. In epithelial cells, required for

产品名称: Ezrin Rabbit Monoclonal Antibody
产品货号: AMRe02993




the formation of microvilli and membrane ruffles on the apical pole. Along with PLEKHG6, required for normal macropinocytosis.


研究领域 (Research Area)

Cell Biology


图片 (Image Data)



Western blot analysis of Ezrin in MEF lysates using Ezrin antibody.



Western blot analysis of Ezrin in HeLa, A549, HL-60, U2OS, C6 lysates using Ezrin antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Ezrin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: Ezrin Rabbit Monoclonal Antibody
产品货号: AMRe02993



注意事项 (Note)

For research use only .