

产品名称: FGFR3 Rabbit Monoclonal Antibody
产品货号: AMRe02998



产品概述 (Summary)

产品名称 (Production Name)	FGFR3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FGFR3
别名 (Alternative Names)	FGFR3; JTK4; Fibroblast growth factor receptor 3; FGFR-3; CD antigen CD333
基因 ID (Gene ID)	2261
蛋白 ID (SwissProt ID)	P22607.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 88 kDa; Observed MW: 125 kDa

研究背景 (Background)

Receptor for acidic and basic fibroblast growth factors. Preferentially binds FGF1.


产品名称: FGFR3 Rabbit Monoclonal Antibody
产品货号: AMRe02998




研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Western blot analysis of FGFR3 in mouse kidney lysates using FGFR3 antibody.



Western blot analysis of FGFR3 in Hela, A549, HL-60, U2OS, C6 lysates using FGFR3 antibody.

注意事项 (Note)

For research use only .