

产品名称: FKBP52 Rabbit Monoclonal Antibody
产品货号: AMRe03000



产品概述 (Summary)

产品名称 (Production Name)	FKBP52 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FKBP4
别名 (Alternative Names)	HBI; p52; Hsp56; FKBP51; FKBP52; FKBP59; PPlase
基因 ID (Gene ID)	2288
蛋白 ID (SwissProt ID)	Q02790.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 52 kDa; Observed MW: 52 kDa

研究背景 (Background)

Immunophilin protein with PPlase and co-chaperone activities (By similarity). Component of unliganded steroid receptors

产品名称: FKBP52 Rabbit Monoclonal Antibody
产品货号: AMRe03000




heterocomplexes through interaction with heat-shock protein 90 (HSP90).


研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Western blot analysis of FKBP52 in Raw264.7 lysates using FKBP52 antibody.



Western blot analysis of FKBP52 in U2OS, C6 lysates using FKBP52 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using FKBP52 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: FKBP52 Rabbit Monoclonal Antibody
产品货号: AMRe03000



注意事项 (Note)

For research use only .