

产品名称: FKBP6 Rabbit Monoclonal Antibody
产品货号: AMRe03001



产品概述 (Summary)

产品名称 (Production Name)	FKBP6 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FKBP6
别名 (Alternative Names)	FKBP36
基因 ID (Gene ID)	8468
蛋白 ID (SwissProt ID)	O75344.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 37 kDa; Observed MW: 37 kDa

研究背景 (Background)

Co-chaperone required during spermatogenesis to repress transposable elements and prevent their mobilization, which is

产品名称: FKBP6 Rabbit Monoclonal Antibody
产品货号: AMRe03001




essential for the germline integrity.


研究领域 (Research Area)

Immunology


图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded mouse testis using FKBP6 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of FKBP6 in rat brain, 3T3, mouse brain, mouse lung lysates using FKBP6 antibody.



Western blot analysis of FKBP6 in HeLa, A549, HL-60, U251, U87-MG lysates using FKBP6 antibody

产品名称: FKBP6 Rabbit Monoclonal Antibody
产品货号: AMRe03001



注意事项 (Note)

For research use only .