

产品名称: FXR1 Rabbit Monoclonal Antibody
产品货号: AMRe03003



产品概述 (Summary)

产品名称 (Production Name)	FXR1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FXR1
别名 (Alternative Names)	FXR1; FXR1P; hFXR1p
基因 ID (Gene ID)	8087
蛋白 ID (SwissProt ID)	P51114.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 70 kDa; Observed MW: 70-80,82-84 kDa

研究背景 (Background)

RNA-binding protein required for embryonic and postnatal development of muscle tissue. May regulate intracellular

产品名称: FXR1 Rabbit Monoclonal Antibody
产品货号: AMRe03003




transport and local translation of certain mRNAs.


研究领域 (Research Area)

Immunology


图片 (Image Data)



Immunocytochemistry analysis of FXR1 (green) in HeLa using FXR1 antibody, and DAPI (blue).




Western blot analysis of FXR1 in PC-12, C2C12 lysates using FXR1 antibody.



Western blot analysis of FXR1 in HeLa, MCF-7, HepG2, U2OS, U251, HL-60 lysates using FXR1 antibody

产品名称: FXR1 Rabbit Monoclonal Antibody
产品货号: AMRe03003



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using FXR1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .