

产品名称: GC1q R Rabbit Monoclonal Antibody
产品货号: AMRe03009



产品概述 (Summary)

产品名称 (Production Name)	GC1q R Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	C1QBP C1QBP; GC1QBP; HABP1; SF2P32; Complement component 1 Q subcomponent-binding protein; mitochondrial; GC1q-R protein; Glycoprotein gC1qBP; C1qBP; Hyaluronan-binding protein 1; Mitochondrial matrix protein p32; p33
别名 (Alternative Names)	
基因 ID (Gene ID)	708
蛋白 ID (SwissProt ID)	Q07021.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 31 kDa; Observed MW: 31 kDa

产品名称: GC1q R Rabbit Monoclonal Antibody
产品货号: AMRe03009




研究背景 (Background)

Binds to the globular "heads" of C1Q thus inhibiting C1 activation.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of GC1q R in A549, HL-60 lysates using GC1q R antibody.



Immunocytochemistry analysis of GC1q R (green) in hela using GC1q R antibody, and DAPI (blue)

注意事项 (Note)

For research use only .