

产品名称: GRSF1 Rabbit Monoclonal Antibody
产品货号: AMRe03016



产品概述 (Summary)

产品名称 (Production Name)	GRSF1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GRSF1
别名 (Alternative Names)	G-rich sequence factor 1
基因 ID (Gene ID)	2926
蛋白 ID (SwissProt ID)	Q12849.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 53 kDa; Observed MW: 53 kDa

研究背景 (Background)

Regulator of post-transcriptional mitochondrial gene expression, required for assembly of the mitochondrial ribosome and

产品名称: GRSF1 Rabbit Monoclonal Antibody
产品货号: AMRe03016




for recruitment of mRNA and lncRNA.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Immunocytochemistry analysis of GRSF1 (green) in CEM using GRSF1 antibody, and DAPI (blue).



Western blot analysis of GRSF1 in HeLa, Ramos lysates using GRSF1 antibody.

注意事项 (Note)

For research use only .