

产品名称: hHR23A Rabbit Monoclonal Antibody
产品货号: AMRe03022



产品概述 (Summary)

产品名称 (Production Name)	hHR23A Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RAD23A
别名 (Alternative Names)	RAD23A; UV excision repair protein RAD23 homolog A; HR23A; hHR23A
基因 ID (Gene ID)	5886
蛋白 ID (SwissProt ID)	P54725.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 40 kDa; Observed MW: 52 kDa

研究背景 (Background)

Multiubiquitin chain receptor involved in modulation of proteasomal degradation. Binds to 'Lys-48'-linked polyubiquitin

产品名称: hHR23A Rabbit Monoclonal Antibody
产品货号: AMRe03022




chains in a length-dependent manner and with a lower affinity to 'Lys-63'-linked polyubiquitin chains.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling; DNA / RNA; DNA Damage & Repair; Base Excision Repair

图片 (Image Data)



Western blot analysis of hHR23A in 3T3 lysates using hHR23A antibody.



Western blot analysis of hHR23A in A549, HL-60, U2OS lysates using hHR23A antibody.

注意事项 (Note)

For research use only .