

产品名称: HIF1 beta Rabbit Monoclonal Antibody
产品货号: AMRe03025



产品概述 (Summary)

产品名称 (Production Name)	HIF1 beta Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat,Hamster

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ARNT ARNT; BHLHE2; Aryl hydrocarbon receptor nuclear translocator; ARNT protein;
别名 (Alternative Names)	Class E basic helix-loop-helix protein 2; bHLHe2; Dioxin receptor; nuclear translocator; Hypoxia-inducible factor 1-beta; HIF-1-beta; HIF1-beta
基因 ID (Gene ID)	405
蛋白 ID (SwissProt ID)	P27540.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 87 kDa; Observed MW: 87 kDa

产品名称: HIF1 beta Rabbit Monoclonal Antibody
产品货号: AMRe03025




研究背景 (Background)

HIF-1 β is also known as AhR nuclear translocator (ARNT) due to its ability to partner with the aryl hydrocarbon receptor (AhR) to form a heterodimeric transcription factor complex. Together with AhR, HIF-1 β plays an important role in xenobiotics metabolism. In addition, a chromosomal translocation leading to a TEL-ARNT fusion protein is associated with acute myeloblastic leukemia. Studies also found that ARNT/HIF-1 β expression levels decrease significantly in pancreatic islets from patients with type 2 diabetes, suggesting that HIF-1 β plays an important role in pancreatic β -cell function.


研究领域 (Research Area)

Cardiovascular

图片 (Image Data)




Western blot analysis of HIF1 beta in Raw264.7, 3T3 lysates using HIF1 beta antibody.



Western blot analysis of HIF1 beta in HeLa, CHO-K1, C6, Ramos lysates using HIF1 beta antibody.

产品名称: HIF1 beta Rabbit Monoclonal Antibody
产品货号: AMRe03025



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using HIF1 beta antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .