

产品名称: HJURP Rabbit Monoclonal Antibody
产品货号: AMRe03030



产品概述 (Summary)

产品名称 (Production Name)	HJURP Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HJURP
别名 (Alternative Names)	FAKTS; URLC9; hFLEG1
基因 ID (Gene ID)	55355
蛋白 ID (SwissProt ID)	Q8NCD3.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 84 kDa; Observed MW: 84 kDa

研究背景 (Background)

Centromeric protein that plays a central role in the incorporation and maintenance of histone H3-like variant CENPA at


产品名称: HJURP Rabbit Monoclonal Antibody
产品货号: AMRe03030




centromeres.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of HJURP in HeLa, CHO-K1 lysates using HJURP antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using HJURP antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .