

产品名称: HLF Rabbit Monoclonal Antibody
产品货号: AMRe03032



产品概述 (Summary)

产品名称 (Production Name)	HLF Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HLF
别名 (Alternative Names)	Hepatic leukemia factor
基因 ID (Gene ID)	3131
蛋白 ID (SwissProt ID)	Q16534.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 33 kDa; Observed MW: 33 kDa

研究背景 (Background)

This gene encodes a member of the proline and acidic-rich (PAR) protein family, a subset of the bZIP transcription factors.


产品名称: HLF Rabbit Monoclonal Antibody
产品货号: AMRe03032



The encoded protein forms homodimers or heterodimers with other PAR family members and binds sequence-specific promoter elements to activate transcription. Chromosomal translocations fusing portions of this gene with the E2A gene cause a subset of childhood B-lineage acute lymphoid leukemias. Alternatively spliced transcript variants have been described, but their biological validity has not been determined.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of HLF in HeLa, A549, HL-60, U87-MG lysates using HLF antibody.

注意事项 (Note)

For research use only .