

产品名称: HMG4 Rabbit Monoclonal Antibody
产品货号: AMRe03033



产品概述 (Summary)

产品名称 (Production Name)	HMG4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HMGB3
别名 (Alternative Names)	HMG4; HMG-4; HMG2A; HMG-2a
基因 ID (Gene ID)	3149
蛋白 ID (SwissProt ID)	O15347.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 23 kDa; Observed MW: 23 kDa

研究背景 (Background)

Multifunctional protein with various roles in different cellular compartments. May act in a redox sensitive manner.

产品名称: HMG4 Rabbit Monoclonal Antibody
产品货号: AMRe03033




Associates with chromatin and binds DNA with a preference to non-canonical DNA structures such as single-stranded DNA. Can bent DNA and enhance DNA flexibility by looping thus providing a mechanism to promote activities on various gene promoters . Proposed to be involved in the innate immune response to nucleic acids by acting as a cytoplasmic promiscuous immunogenic DNA/RNA sensor . Negatively regulates B-cell and myeloid cell differentiation. In hematopoietic stem cells may regulate the balance between self-renewal and differentiation. Involved in negative regulation of canonical Wnt signaling .


研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)




Western blot analysis of HMG4 in 3T3 lysates using HMG4 antibody.



Western blot analysis of HMG4 in HeLa, A549, HL-60, U2OS, C6 lysates using HMG4 antibody.

产品名称: HMG4 Rabbit Monoclonal Antibody
产品货号: AMRe03033



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using HMG4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .