

产品名称: HspA14 Rabbit Monoclonal Antibody
产品货号: AMRe03038



产品概述 (Summary)

产品名称 (Production Name)	HspA14 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HSPA14
别名 (Alternative Names)	HSPA14; HSP60; HSP70L1; Heat shock 70 kDa protein 14; HSP70-like protein 1; Heat shock protein HSP60
基因 ID (Gene ID)	51182
蛋白 ID (SwissProt ID)	Q0VDF9.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 55 kDa; Observed MW: 55 kDa

研究背景 (Background)

产品名称: HspA14 Rabbit Monoclonal Antibody
产品货号: AMRe03038




Component of the ribosome-associated complex (RAC), a complex involved in folding or maintaining nascent polypeptides in a folding-competent state. In the RAC complex, binds to the nascent polypeptide chain, while DNAJC2 stimulates its ATPase activity.


研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Immunocytochemistry analysis of HspA14 (green) in Jurkat using HspA14 antibody, and DAPI (blue).



Western blot analysis of HSPA14 in 3T3 lysates using HSPA14 antibody.



Western blot analysis of HSPA14 in HL-60, U251, U87-MG lysates using HSPA14 antibody

产品名称: HspA14 Rabbit Monoclonal Antibody
产品货号: AMRe03038



注意事项 (Note)

For research use only .