

产品名称: ING1 Rabbit Monoclonal Antibody
产品货号: AMRe03049



产品概述 (Summary)

产品名称 (Production Name)	ING1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ING1
别名 (Alternative Names)	p33; p47; p33ING1; p24ING1c; p33ING1b; p47ING1a
基因 ID (Gene ID)	3621
蛋白 ID (SwissProt ID)	Q9UK53.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 47 kDa; Observed MW: 47 kDa

研究背景 (Background)

This gene encodes a tumor suppressor protein that can induce cell growth arrest and apoptosis. The encoded protein is a

产品名称: ING1 Rabbit Monoclonal Antibody
产品货号: AMRe03049




nuclear protein that physically interacts with the tumor suppressor protein TP53 and is a component of the p53 signaling pathway. Reduced expression and rearrangement of this gene have been detected in various cancers. Multiple alternatively spliced transcript variants encoding distinct isoforms have been reported.


研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of ING1 in C6, Raw264.7 lysates using ING1 antibody.



Western blot analysis of ING1 in Hela, A549, HL-60, U251, U87-MG lysates using ING1 antibody.

注意事项 (Note)

For research use only .