

产品名称: KIFAP3 Rabbit Monoclonal Antibody
产品货号: AMRe03057



产品概述 (Summary)

产品名称 (Production Name)	KIFAP3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KIFAP3
别名 (Alternative Names)	KAP3; KIF3AP; Kifap3; Kinesin associated protein 3; SMAP
基因 ID (Gene ID)	22920
蛋白 ID (SwissProt ID)	Q92845.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 91 kDa; Observed MW: 91 kDa

研究背景 (Background)

Involved in tethering the chromosomes to the spindle pole and in chromosome movement. Binds to the tail domain of the

产品名称: KIFAP3 Rabbit Monoclonal Antibody
产品货号: AMRe03057




KIF3A/KIF3B heterodimer to form a heterotrimeric KIF3 complex and may regulate the membrane binding of this complex.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of KIFAP3 in A549, HL-60, C6 lysates using KIFAP3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using KIFAP3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .