

**产品名称: Ku70 Rabbit Monoclonal Antibody**  
**产品货号: AMRe03060**



## 产品概述 (Summary)

|                        |                                        |
|------------------------|----------------------------------------|
| 产品名称 (Production Name) | Ku70 Rabbit Monoclonal Antibody        |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody |
| 宿主 (Host)              | Rabbit                                 |
| 应用 (Application)       | WB,IHC,IP                              |
| 种属反应性 (Reactivity)     | Human                                  |

## 产品性能 (Performance)

|                     |                                                                                                      |
|---------------------|------------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                         |
| 修饰 (Modification)   | Unmodified                                                                                           |
| 同种型 (Isotype)       | IgG                                                                                                  |
| 克隆 (Clonality)      | Monoclonal                                                                                           |
| 形式 (Form)           | Liquid                                                                                               |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.             |
| 储存溶液 (Buffer)       | 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein |
| 纯化方式 (Purification) | Affinity Purification                                                                                |

## 免疫原信息 (Immunogen)

|                        |                                                                                                                                                                                                                                                           |
|------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | XRCC6<br>XRCC6; G22P1; X-ray repair cross-complementing protein 6; 5'-deoxyribose-5-phosphate lyase Ku70; 5'-dRP lyase Ku70; 70 kDa subunit of Ku antigen; ATP-dependent DNA helicase 2 subunit 1; ATP-dependent DNA helicase II 70 kDa subunit; CTC box- |
| 别名 (Alternative Names) |                                                                                                                                                                                                                                                           |
| 基因 ID (Gene ID)        | 2547                                                                                                                                                                                                                                                      |
| 蛋白 ID (SwissProt ID)   | P12956.                                                                                                                                                                                                                                                   |

## 产品应用 (Application)

|                          |                                             |
|--------------------------|---------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50 |
| 蛋白分子量 (Molecular Weight) | Calculated MW: 70 kDa; Observed MW: 70 kDa  |

产品名称: Ku70 Rabbit Monoclonal Antibody  
产品货号: AMRe03060




## 研究背景 (Background)

It works in the 3'-5' direction. Binding to DNA may be mediated by XRCC6. Involved in DNA non-homologous end joining (NHEJ) required for double-strand break repair and V(D)J recombination. The XRCC5/6 dimer acts as regulatory subunit of the DNA-dependent protein kinase complex DNA-PK by increasing the affinity of the catalytic subunit PRKDC to DNA by 100-fold. The XRCC5/6 dimer is probably involved in stabilizing broken DNA ends and bringing them together.


## 研究领域 (Research Area)

Epigenetics and Nuclear Signaling

## 图片 (Image Data)



Western blot analysis of Ku70 in HeLa, Ramos lysates using Ku70 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Ku70 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

## 注意事项 (Note)

For research use only .