

产品名称: LIG1 Rabbit Monoclonal Antibody
产品货号: AMRe03065



产品概述 (Summary)

产品名称 (Production Name)	LIG1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LIG1
别名 (Alternative Names)	LIG1; DNA ligase 1; DNA ligase I; Polydeoxyribonucleotide synthase [ATP] 1
基因 ID (Gene ID)	3978
蛋白 ID (SwissProt ID)	P18858.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 102 kDa; Observed MW: 102 kDa

研究背景 (Background)

DNA ligase that seals nicks in double-stranded DNA during DNA replication, DNA recombination and DNA repair.


产品名称: **LIG1 Rabbit Monoclonal Antibody**
产品货号: **AMRe03065**




研究领域 (Research Area)

Epigenetics and Nuclear Signaling


图片 (Image Data)



Immunocytochemistry analysis of LIG1 (green) in HeLa using LIG1 antibody, and DAPI (blue).




Western blot analysis of LIG1 in PC-12 lysates using LIG1 antibody.



Western blot analysis of LIG1 in HeLa, A549, HL-60 lysates using LIG1 antibody

产品名称: LIG1 Rabbit Monoclonal Antibody
产品货号: AMRe03065



Immunohistochemistry analysis of paraffin-embedded Human tonsil using LIG1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .