

产品名称: MLLT6 Rabbit Monoclonal Antibody
产品货号: AMRe03090



产品概述 (Summary)

产品名称 (Production Name)	MLLT6 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MLLT6
别名 (Alternative Names)	AF17
基因 ID (Gene ID)	4302
蛋白 ID (SwissProt ID)	P55198.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 112 kDa; Observed MW: 112 kDa

研究背景 (Background)

Enables histone binding activity and nucleosome binding activity. Predicted to be involved in positive regulation of


产品名称: MLLT6 Rabbit Monoclonal Antibody
产品货号: AMRe03090



transcription by RNA polymerase II. Predicted to act upstream of or within several processes, including negative regulation of histone H3-K79 methylation; renal potassium excretion; and renal sodium excretion. Predicted to be active in nucleus.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of MLLT6 in HeLa, A549, HL-60, U2OS lysates using MLLT6 antibody.

注意事项 (Note)

For research use only .