

产品名称: Interferon Regulatory Factor 4 Rabbit Monoclonal Antibody
产品货号: AMRe03098



产品概述 (Summary)

产品名称 (Production Name)	Interferon Regulatory Factor 4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IRF4
别名 (Alternative Names)	IRF4; MUM1; Interferon regulatory factor 4; IRF-4; Lymphocyte-specific interferon regulatory factor; LSIRF; Multiple myeloma oncogene 1; NF-EM5
基因 ID (Gene ID)	3662
蛋白 ID (SwissProt ID)	Q15306.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 52 kDa; Observed MW: 52 kDa

产品名称: Interferon Regulatory Factor 4 Rabbit Monoclonal Antibody
产品货号: AMRe03098




研究背景 (Background)

Melanoma associated antigen (mutated) 1 (MUM1, EXPAND1) is a PWWP-domain containing chromatin binding protein involved in maintaining chromatin architecture of interphase chromosomes.

研究领域 (Research Area)

Tags & Cell Markers

图片 (Image Data)



Western blot analysis of MUM1 in Ramos lysates using Interferon Regulatory Factor 4 antibody.

注意事项 (Note)

For research use only .