

产品名称: NME2 Rabbit Monoclonal Antibody
产品货号: AMRe03114



产品概述 (Summary)

产品名称 (Production Name)	NME2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NME2 NME2; NM23B; Nucleoside diphosphate kinase B; NDK B; NDP kinase B; C-myc
别名 (Alternative Names)	purine-binding transcription factor PUF; Histidine protein kinase NDKB; nm23-H2
基因 ID (Gene ID)	4831
蛋白 ID (SwissProt ID)	P22392.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 17 kDa; Observed MW: 17 kDa

产品名称: NME2 Rabbit Monoclonal Antibody
产品货号: AMRe03114




研究背景 (Background)

Major role in the synthesis of nucleoside triphosphates other than ATP. Negatively regulates Rho activity by interacting with AKAP13/LBC.


研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Western blot analysis of NME2 in mouse spleen lysates using NME2 antibody.



Western blot analysis of NME2 in U2OS, C6, Ramos lysates using NME2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using NME2 antibody. High-pressure and

产品名称: NME2 Rabbit Monoclonal Antibody
产品货号: AMRe03114



temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .