

产品名称: MYH10 Rabbit Monoclonal Antibody
产品货号: AMRe03117



产品概述 (Summary)

产品名称 (Production Name)	MYH10 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MYH10
别名 (Alternative Names)	NMMHCB; NMMHC-IIB
基因 ID (Gene ID)	4628
蛋白 ID (SwissProt ID)	P35580.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 229 kDa; Observed MW: 229 kDa

研究背景 (Background)

Cellular myosin that appears to play a role in cytokinesis, cell shape, and specialized functions such as secretion and

产品名称: MYH10 Rabbit Monoclonal Antibody
产品货号: AMRe03117




capping.


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of nonmuscle Myosin IIB/MYH10 in 3T3 lysates using MYH10 antibody.



Western blot analysis of nonmuscle Myosin IIB/MYH10 in A549, HL-60, U251, U87-MG lysates using nonmuscle Myosin IIB/MYH10 antibody.

注意事项 (Note)

For research use only .