

产品名称: NSDHL Rabbit Monoclonal Antibody
产品货号: AMRe03120



产品概述 (Summary)

产品名称 (Production Name)	NSDHL Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NSDHL
别名 (Alternative Names)	H105E3; XAP104; SDR31E1
基因 ID (Gene ID)	50814
蛋白 ID (SwissProt ID)	Q15738.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 42 kDa; Observed MW: 42 kDa

研究背景 (Background)

The protein encoded by this gene is localized in the endoplasmic reticulum and is involved in cholesterol biosynthesis.

产品名称: NSDHL Rabbit Monoclonal Antibody
产品货号: AMRe03120




Mutations in this gene are associated with CHILD syndrome, which is a X-linked dominant disorder of lipid metabolism with disturbed cholesterol biosynthesis, and typically lethal in males. Alternatively spliced transcript variants with differing 5' and 3' UTR have been found for this gene.

研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Western blot analysis of NSDHL in HeLa, A549, HL-60, U251 lysates using NSDHL antibody.

注意事项 (Note)

For research use only .