

产品名称: NUP107 Rabbit Monoclonal Antibody
产品货号: AMRe03122



产品概述 (Summary)

产品名称 (Production Name)	NUP107 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NUP107
别名 (Alternative Names)	ODG6; NUP84; NPHS11; ODG6; GAMOS7
基因 ID (Gene ID)	57122
蛋白 ID (SwissProt ID)	P57740.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 106 kDa; Observed MW: 106 kDa

研究背景 (Background)

Plays a role in the nuclear pore complex (NPC) assembly and/or maintenance (PubMed:12552102, PubMed:15229283).

产品名称: NUP107 Rabbit Monoclonal Antibody
产品货号: AMRe03122




Required for the assembly of peripheral proteins into the NPC (PubMed:15229283, PubMed:12552102). May anchor NUP62 to the NPC (PubMed:15229283).

研究领域 (Research Area)

Tags & Cell Markers

图片 (Image Data)



Western blot analysis of Nup107 in A549, HL-60, U251 lysates using NUP107 antibody.

注意事项 (Note)

For research use only .