

产品名称: NUP50 Rabbit Monoclonal Antibody
产品货号: AMRe03125



产品概述 (Summary)

产品名称 (Production Name)	NUP50 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NUP50
别名 (Alternative Names)	NPAP60; NPAP60L
基因 ID (Gene ID)	10762
蛋白 ID (SwissProt ID)	Q9UKX7.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 50 kDa; Observed MW: 50 kDa

研究背景 (Background)


Component of the nuclear pore complex that has a direct role in nuclear protein import.

产品名称: NUP50 Rabbit Monoclonal Antibody
产品货号: AMRe03125




研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of NUP50 in HeLa, C6 lysates using NUP50 antibody.



Immunocytochemistry analysis of NUP50 (green) in HeLa using NUP50 antibody, and DAPI (blue)

注意事项 (Note)

For research use only.