

产品名称: ORAI3 Rabbit Monoclonal Antibody
产品货号: AMRe03129



产品概述 (Summary)

产品名称 (Production Name)	ORAI3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ORAI3
别名 (Alternative Names)	TMEM142C
基因 ID (Gene ID)	93129
蛋白 ID (SwissProt ID)	Q9BRQ5.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 31 kDa; Observed MW: 31 kDa

研究背景 (Background)

Ca²⁺ release-activated Ca²⁺-like (CRAC-like) channel subunit which mediates Ca²⁺ influx and increase in Ca²⁺-selective


产品名称: ORAI3 Rabbit Monoclonal Antibody
产品货号: AMRe03129




current by synergy with the Ca²⁺ sensor, STIM1.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Orai3 in HeLa, A549, HL-60, U251, U87-MG lysates using ORAI3 antibody.



Immunohistochemistry analysis of paraffin-embedded liver cancer using Orai3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .