

产品名称: p27 Kip 1 Rabbit Monoclonal Antibody
产品货号: AMRe03132



产品概述 (Summary)

产品名称 (Production Name)	p27 Kip 1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Cdkn1B
别名 (Alternative Names)	p27; Kip1; p27Kip1; AA408329; AI843786
基因 ID (Gene ID)	12576.0
蛋白 ID (SwissProt ID)	P46414.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 22 kDa; Observed MW: 27 kDa

研究背景 (Background)

Important regulator of cell cycle progression. Inhibits the kinase activity of CDK2 bound to cyclin A, but has little inhibitory

产品名称: p27 Kip 1 Rabbit Monoclonal Antibody
产品货号: AMRe03132




activity on CDK2 bound to SPDYA (By similarity).


研究领域 (Research Area)

Cell Biology


图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded mouse lung using p27 Kip 1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of p27 KIP 1 in 3T3 lysates using p27 KIP 1 antibody.



Western blot analysis of p27 KIP 1 in HeLa, A549, HL-60 lysates using p27 KIP 1 antibody

产品名称: p27 Kip 1 Rabbit Monoclonal Antibody
产品货号: AMRe03132



注意事项 (Note)

For research use only .