

产品名称: P2Y13 Rabbit Monoclonal Antibody
产品货号: AMRe03134



产品概述 (Summary)

产品名称 (Production Name)	P2Y13 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	P2RY13
别名 (Alternative Names)	GPCR1; GPR86; GPR94; P2Y13; SP174; FKSG77
基因 ID (Gene ID)	53829
蛋白 ID (SwissProt ID)	Q9BPV8.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 41 kDa; Observed MW: 41 kDa

研究背景 (Background)

The product of this gene belongs to the family of G-protein coupled receptors. This family has several receptor subtypes

产品名称: P2Y13 Rabbit Monoclonal Antibody
产品货号: AMRe03134




with different pharmacological selectivity, which overlaps in some cases, for various adenosine and uridine nucleotides. This receptor is activated by ADP. [provided by RefSeq, Sep 2008]

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of P2Y13 in HeLa, A549 lysates using P2Y13 antibody.

注意事项 (Note)

For research use only .