

产品名称: PELP1 Rabbit Monoclonal Antibody
产品货号: AMRe03139



产品概述 (Summary)

产品名称 (Production Name)	PELP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PELP1
别名 (Alternative Names)	HMX3; MNAR; P160; PELP1; PELP1 proline glutamic acid leucine rich protein 1; PELP1 proline- glutamic
基因 ID (Gene ID)	27043
蛋白 ID (SwissProt ID)	Q8IZL8.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 120 kDa; Observed MW: 160 kDa

研究背景 (Background)

产品名称: PELP1 Rabbit Monoclonal Antibody
产品货号: AMRe03139




This gene encodes a transcription factor which coactivates transcription of estrogen receptor responsive genes and corepresses genes activated by other hormone receptors or sequence-specific transcription factors. Expression of this gene is regulated by both members of the estrogen receptor family. This gene may be involved in the progression of several types of cancer. Alternative splicing results in multiple transcript variants.


研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Immunocytochemistry analysis of PELP1 (green) in HeLa using PELP1 antibody, and DAPI (blue).



Western blot analysis of PELP1 in HeLa, MCF-7, HepG2 lysates using PELP1 antibody.



产品名称: PELP1 Rabbit Monoclonal Antibody
产品货号: AMRe03139



Immunohistochemistry analysis of paraffin-embedded Human breast cancer tissue using PELP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .