

产品名称: Poliovirus Receptor Rabbit Monoclonal Antibody
产品货号: AMRe03156



产品概述 (Summary)

产品名称 (Production Name)	Poliovirus Receptor Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PVR
别名 (Alternative Names)	PVR; PVS; Poliovirus receptor; Nectin-like protein 5; NECL-5; CD155
基因 ID (Gene ID)	5817
蛋白 ID (SwissProt ID)	P15151.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 45 kDa; Observed MW: 60-80 kDa

研究背景 (Background)

Mediates NK cell adhesion and triggers NK cell effector functions. Binds two different NK cell receptors: CD96 and CD226.

产品名称: Poliovirus Receptor Rabbit Monoclonal Antibody
产品货号: AMRe03156




These interactions accumulates at the cell-cell contact site, leading to the formation of a mature immunological synapse between NK cell and target cell.

研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of Poliovirus Receptor/PVR in HeLa, A549, U2OS lysates using Poliovirus Receptor antibody.

注意事项 (Note)

For research use only .