

产品名称: Protein CASP Rabbit Monoclonal Antibody
产品货号: AMRe03165



产品概述 (Summary)

产品名称 (Production Name)	Protein CASP Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CUX1
别名 (Alternative Names)	CDP; CUX; p75; CASP; CDP1; COY1; Clox; GDDI; p100; p110; p200; CUTL1; GOLIM6; CDP/Cut; Cux/CDP; Nbla10317
基因 ID (Gene ID)	1523
蛋白 ID (SwissProt ID)	Q13948.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 77 kDa; Observed MW: 77 kDa

研究背景 (Background)


产品名称: Protein CASP Rabbit Monoclonal Antibody
产品货号: AMRe03165




May be involved in intra-Golgi retrograde transport.

研究领域 (Research Area)


图片 (Image Data)



Immunocytochemistry analysis of Protein CASP (green) in A549 using Protein CASP antibody, and DAPI (blue).



Western blot analysis of Protein CASP in 3T3, Raw264.7 lysates using Protein CASP antibody.



Western blot analysis of Protein CASP in HeLa, A549, HL-60, U2OS, C6 lysates using Protein CASP antibody

产品名称: Protein CASP Rabbit Monoclonal Antibody
产品货号: AMRe03165



注意事项 (Note)

For research use only .