

产品名称: PRPF4 Rabbit Monoclonal Antibody
产品货号: AMRe03166



产品概述 (Summary)

产品名称 (Production Name)	PRPF4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PRPF4
别名 (Alternative Names)	PRP4; RP70; HPRP4; Prp4p; HPRP4P; SNRNP60
基因 ID (Gene ID)	9128
蛋白 ID (SwissProt ID)	O43172.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 58 kDa; Observed MW: 58 kDa

研究背景 (Background)

Participates in pre-mRNA splicing. Part of the U4/U5/U6 tri-snRNP complex, one of the building blocks of the spliceosome.


产品名称: PRPF4 Rabbit Monoclonal Antibody
产品货号: AMRe03166




研究领域 (Research Area)

Cell Biology


图片 (Image Data)



Western blot analysis of PRPF4 in Raw264.7, C6 lysates using PRPF4 antibody.



Western blot analysis of PRPF4 in HeLa, A549, HL-60 lysates using PRPF4 antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using PRPF4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: PRPF4 Rabbit Monoclonal Antibody
产品货号: AMRe03166



注意事项 (Note)

For research use only .