

产品名称: PTBP1 Rabbit Monoclonal Antibody
产品货号: AMRe03168



产品概述 (Summary)

产品名称 (Production Name)	PTBP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Hamster

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PTBP1
别名 (Alternative Names)	PTB; PTB2; PTB3; PTB4; pPTB; HNRPI; PTB-1; PTB-T; HNRNPI; HNRNP-I
基因 ID (Gene ID)	5725
蛋白 ID (SwissProt ID)	P26599.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 57 kDa; Observed MW: 57 kDa

研究背景 (Background)

Plays a role in pre-mRNA splicing and in the regulation of alternative splicing events. Binds to the polypyrimidine tract of

产品名称: PTBP1 Rabbit Monoclonal Antibody
产品货号: AMRe03168




introns. May promote RNA looping when bound to two separate polypyrimidine tracts in the same pre-mRNA.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling


图片 (Image Data)



Western blot analysis of PTBP1 in 3T3 lysates using PTBP1 antibody.



Western blot analysis of PTBP1 in HeLa, CHO-K1, Ramos lysates using PTBP1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using PTBP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: PTBP1 Rabbit Monoclonal Antibody
产品货号: AMRe03168



注意事项 (Note)

For research use only .